

# Raman Analysis to Evaluate Carbon Coating on Graphite Anodes in LIBs

## Presenting Author

Hajar Elazri

HORIBA FRANCE SAS, 14 boulevard Thomas Gobert, Palaiseau, France

[hajar.elazri@horiba.com](mailto:hajar.elazri@horiba.com)

Enping Hu

HORIBA(CHINA)TRADING CO.,LTD, 1068 West Tianshan Road, Shanghai, China


[enping.hu@horiba.com](mailto:enping.hu@horiba.com)

## Abstract


In lithium-ion battery (LIB) anode materials, graphite is commonly coated with amorphous carbon to enhance lithium capacity. Raman spectroscopy is a powerful tool for assessing coating thickness and uniformity, as shown in Figure 1. A thicker coating corresponds to a higher  $I_D/I_G$  ratio. As particles coating are usually non-uniform, statistical analysis of hundreds of particles is required to evaluate the overall coating quality.

HORIBA provides two methods for this statistical analysis. The first approach utilizes ParticleFinder with a vacuum facility to disperse the powder into individual particles, enabling Raman measurements of each particle and determining the mean  $I_D/I_G$  value. The second approach employs Raman mapping to visualize  $I_D/I_G$  distribution. With the help of HORIBA QCarbon, the ratio distribution, its mean value, and standard deviation could be obtained by one click, streamlining the analysis process.

## Figures



**Figure 1:** Left: Raman spectra of **graphite**, **carbon** and **carbon coated graphite** with different thickness; Right: QCarbon app to evaluate  $I_D/I_G$  by one click.



**Figure 2:** Two methods to evaluate the coating degree. Upper: with ParticleFinder method; lower: with mapping method.