3D sensing platform for single cell extracellular pH mapping in time and space: a pre-clinical model to study tumor microenvironment and drug screening

Stefania Forciniti1,

Valentina Onesto<sup>1</sup>, Riccardo Rizzo<sup>1</sup>, Anil Chandra<sup>1</sup> Saumya Prasad<sup>1</sup>, Helena Iuele<sup>1</sup>, Francesco Colella<sup>1</sup>, Giuseppe Gigli<sup>1,2</sup>, Loretta L. del Mercato<sup>1</sup>

<sup>1</sup>Institute of Nanotechnology, National Research Council (CNR-NANOTEC), c/o Campus Ecotekne, via Monteroni, 73100, Lecce, Italy;

<sup>2</sup>Department of Mathematics and Physics "Ennio De Giorgi", University of Salento, via Arnesano, 73100, Lecce, Italy.

Stefania.forciniti@nanotec.cnr.it

## Abstract

Elevated hydrogen ion concentration accumulates in the tumor microenvironment as a result of metabolic reprogramming of cancer cells towards increased aerobic glycolysis and reduced mitochondrial oxidative phosphorylation, referred to as Warburg effect [1]. This results in spatially and temporally heterogenous acidic microenvironment which affects cancer initiation and progression as well as the efficacy of anti-cancer drug treatments [2, 3].

Therefore, monitoring the local pH metabolic fluctuations is crucial in understanding the basic biology of the tumors, and can also be used as a valid metabolic readout for cancer diagnosis and treatment [4].

Here, a method to embed silica-based ratiometric fluorescent pH sensors into a spherical 3D cell culture system, coupled with a computational method for pH spatio-temporal mapping, is presented (see figure 1) [5]. By using a confocal laser fluorescence microscopy, 3D time-lapse imaging of living cells was performed and the extracellular pH variations were monitored between 3D scaffolds with either mono or co-cultures of tumor and stroma pancreatic cells.

We obtained that the extracellular pH is cell linespecific and time-dependent. Indeed, differences in pH were detected between 3D monocultures compared to co-cultures, thus suggesting a metabolic crosstalk between tumor and stroma cells.

In conclusion, the 3D cell culture system has the potential of imaging complex 3D co-coltures in real time and of detecting their pH metabolic interplay under controlled experimental conditions, making it a suitable platform for drug screening and personalized medicine.

The research leading to these results received funding from the European Research Council (ERC) under the European Union's Horizon 2020 research and innovation programme (grant agreement No 759959, ERC-StG "INTERCELLMED").

## References

- Schwartz L, Supuran CT, Alfarouk KO. The Warburg Effect and the Hallmarks of Cancer. Anticancer Agents Med Chem. 2017;17(2):164-170;
- [2] Boedtkjer E, Pedersen SF. The Acidic Tumor Microenvironment as a Driver of Cancer. Annu Rev Physiol. 2020 Feb 10;82:103-126;
- [3] Zhang A, Miao K, Sun H, Deng CX. Tumor heterogeneity reshapes the tumor microenvironment to influence drug resistance. Int J Biol Sci. 2022 Apr 24;18(7):3019-3033;
- Parks SK, Chiche J, Pouyssegur J. pH control mechanisms of tumor survival and growth. J Cell Physiol. 2011 Feb;226(2):299-308;
- [5] Rizzo R, Onesto V, Forciniti S, Chandra A, Prasad S, luele H, Colella F, Gigli G, Del Mercato LL. A pH-sensor scaffold for mapping spatiotemporal gradients in three-dimensional in vitro tumour models. Biosens Bioelectron. 2022 Sep 15;212:114401.

## **Figures**



Figure 1. Schematic illustration of the method to embed silica-based fluorescent pH sensors into alginate 3D microgels tumor models. i) Tumor/stromal cells and sensors microencapsulation in alginate to form 3D microgels tumor models; ii) pH-sensing hydrogel-based tumor/stroma model; iii) Spatio-Temporal extracellular pH mapping obtained with a specific computational analysis.