

MAY 06, 2020
CONFERENCE
ONLINE



ANTIBODY-GATED INDICATOR RELEASING MESOPOROUS MATERIALS: A POTENTIAL BIOSENSOR PLATFORM TO BE USED IN THE DEVELOPMENT OF RAPID TESTS

Estela Climent, Elena Costa, Kornelia Gawlitza, Wei Wan, Michael G. Weller, Knut Rurack


Chemical and Optical Sensing Division, Bundesanstalt für Materialforschung und -prüfung (BAM)

Motivation

Keeping in mind the high sensitivity offered by gated indicator-releasing micro- and nanoparticles due to their inherent features of signal amplification, we embarked on system optimisation to develop a potential biosensor platform for use in rapid tests. Here, the insecticide permethrin, a type-I pyrethroid, was selected as target analyte, because type-I pyrethroids play an important role in airplane disinfection.



Design and synthesis routes of materials



Response in solution

1.- Analytical assay optimisation with S1.1-AB (FLU as dye)




Fig. 1 Increase in fluorescence at 522 nm ($\lambda_{\text{exc}} = 490 \text{ nm}$) in the presence (red) and the absence (black) of permethrin (6.6 mg kg^{-1}) of the supernatants of S1.1-AB suspensions.

2.- Evaluation of pore size and different loading routes followed using SRG as dye






Fig. 2 Increase in fluorescence at 550 nm ($\lambda_{\text{exc}} = 532 \text{ nm}$) in the presence (red) and the absence (black) of 3-PBA (0.5 mg kg^{-1}) of the supernatants of S1.1-AB suspensions.

Materials characterisation

Transmission electron microscopy (TEM)



N₂ Adsorption/desorption isotherms and pore size distribution




Elemental (EA) and thermogravimetric (TGA) analysis


	EA					TGA
	FLU %	SRG %	HAPten %	PEG %	TOTAL %	TOTAL %
S1.1	6.6	—	2.0	—	8.6	9.4
S1.2	—	29.4	2.0	—	31.3	29.6
S2.2a	—	15.0	1.7	15.2	31.9	30.8
S2.2b	—	17.0	3.2	14.8	35.0	34.3
S3.2a	—	12.3	1.6	17.6	30.2	30.0
S4.2a	—	21.4	2.2	11.1	34.7	31.3
S5.2a	—	16.7	3.5	10.2	30.4	29.3

3.- Concentration-dependent studies (after 5 min of reaction)

Antibody at the surface (with different dyes)




Antibody inside of the pores (S2.2a-AB) Similar sensitivity



Signal amplification factor

Solid	C _{10%} μg Kg ⁻¹	AF
S1.1-AB	13	2.5
S1.2-AB	1	860
S2.2a-AB	2.5	1010
S2.2b-AB	4	197
S3.2-AB	4	12
S4.2-AB	—	—
S5.2-AB	—	—

Incorporation into different sensing platforms



Conclusions

- ✓ Increase in ionic strength and presence of organic solvents accelerates release kinetics
- ✓ Indicator dye is decisive for the sensitivity
- ✓ Better material response is observed when size of Fab region of the antibody is similar to the pore diameter and hapten is grafted at the entrance of the pores

CONTACT PERSON

Estela Climent Terol

Chemical and Optical Sensing Division
Bundesanstalt für Materialforschung und -prüfung (BAM)

Estela.climent@bam.de

REFERENCES

- E. Climent, A. Bernardos, R. Martinez-Manez, A. Maquieira, M. D. Marcos, N. Pastor-Navarro, R. Puchades, F. Sancenon, J. Soto, P. Amoros, J. Am. Chem. Soc 131 (2009) 14075.
- J. Bell, E. Climent, M. Hecht, M. Buurman, K. Rurack, ACS Sens. 1 (2016) 334.
- E. Climent, D. Groninger, M. Hecht, M. A. Walter, R. Martinez-Manez, M. G. Weller, F. Sancenon, P. Amoros, K. Rurack, Chem. Eur. J. 19 (2013) 4117.
- E. Costa, E. Climent, S. Ast, M. G. Weller, J. Canning, K. Rurack, Analyst (2020) DOI: 10.1039/D0AN00319K.
- E. Costa, E. Climent, K. Gawlitza, W. Wan, M. G. Weller, K. Rurack, J. Mat. Chem. B (2020) DOI: 10.1039/d0tb00371a.

