

Unravelling the interaction of MXenes with cell membrane and extracellular vesicles

Maha Alhosani^{1*}, Olla Elkhatab^{1*}, Roberta Cagliani^{1*}, Benjamin Chacon², Jawaher Alhammedi³, Laura Fusco^{2,4}, Linda Giro⁴, Marco Orecchioni⁵, Yury Gogotsi², Anna-Maria Pappa^{3,6*}, and Lucia Gemma Delogu^{1,4*}

*Those authors contributed equally to this work

¹Department of Biological Sciences, Khalifa University of Science and Technology, Abu Dhabi, UAE.

²A. J. Drexel Nanomaterials Institute and Department of Materials Science and Engineering, Drexel University, Philadelphia, PA, United States.

³Department of Biomedical Engineering, Khalifa University of Science and Technology, Abu Dhabi, UAE.

⁴Department of Biomedical Science, University of Padua, Padua, Italy.

⁵Immunology Center of Georgia, Augusta University, Augusta, Georgia.

⁶Center of Biotechnology, Khalifa University of Science and Technology, Abu Dhabi, UAE.

100045421@ku.ac.ae

Transition metal carbides, nitrides, and carbonitrides (MXene) are emerging 2D materials with versatile chemistry, garnering significant interest in diagnostics and biomedicine [1]. Recent studies have revealed safe interactions of MXene at tissue and cellular level [2], signalling the potential possibilities for biomedical applications, such as encapsulation within extracellular vesicles (EVs) for future therapeutics. EVs are nanosized, membrane-bound particles that transfer biomolecules such as proteins, lipids, RNAs, and DNA into target cells [3]. Accumulating evidence shows that EVs have a role in intracellular communication, influencing processes like immune response and cancer progression [3,4]. Understanding the chemistry and the intracellular communication between MXene and EVs may open new avenues to revolutionize targeted drug delivery. In this study, MXene-encapsulated EVs are extracted from immune cells of THP-1 cells by mechanical extrusion protocol. The obtained changes in hydrodynamic diameter and surface charges shed light on the dynamics of the interaction of MXene with EVs at different incubation time points. Additionally, flow cytometric analysis is performed to gain deeper insight for the presence of specific extracellular and intracellular EVs surface markers. Future studies aim to delve into the potential use of MXene encapsulated EVs, extracted from THP-1 immune cells, for transfection of biomolecules and therapeutic tagging for biomedical purposes.

References

- [1] Chen, L. et al. Biomedical applications of mxenes: From nanomedicine to Biomaterials. *Accounts of Materials Research* (2022).
- [2] Fusco, L. et al. Immune profiling and multiplexed Label-Free detection of 2D MXenes by mass cytometry and High-Dimensional imaging. *Advanced Materials* (2022).
- [3] Xu, R. et al. Extracellular vesicle isolation and characterization: toward clinical application. *Journal of Clinical Investigation* (2016).
- [4] Balaj, L. et al. Tumour microvesicles contain retrotransposon elements and amplified oncogene sequences. *Nature Communications*, (2011).